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Operating Instructions

Safety Cabinets

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Brandveiligheidsopslagkasten

Mode d'emploi

Armoires anti-feu

Manual de instrucciones

Armarios ignífugos

Deutsch

English

Nederlands

Français

Español

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Contents of these Operating Instructions

1. Models
2. Technical data / Drawings
3. Intended uses
4. General safety information / Warranty
5. Transport
6. Installation
7. Interior fittings
8. Storage amounts / Storing substances together
9. Ventilation
10. Grounding
11. Markings
12. Decommissioning after a fire
13. Care and maintenance
14. Disposal

1. Models

Type 90 – safety cabinets as per EN 14470-1

Type 30 – safety cabinets as per EN 14470-1

G 30 – gas cylinder safety cabinets as per EN 14470-2

2. Technical data / Drawings

see Annex 1

3. Intended uses

- **Type 30 / Type 90 safety cabinets are suitable for the storage of flammable liquids in closed containers at the workplace while observing the relevant national rules and regulations.**
- **G 30 gas cylinder safety cabinets are suitable for the storage of gas cylinders at the workplace while observing the relevant national rules and regulations.**



4. General safety information / Warranty

- For the operation of your safety cabinet you must observe the national rules, regulations and laws that are applicable.
- Before starting to use your safety cabinet you must carry out a functional check in order to recognise any damage that may have been caused by transport, importation or improper installation!
- Safety devices must be checked at regular intervals – this should only be carried out by our trained specialists! This guarantees proper and safe long-term storage.
- Warranty: in accordance with our General Terms and Conditions. See also item 13.
- Please also observe your duty of exercising due care – any damage must be remedied and spilt liquids removed immediately!
- The area in which the doors open must be kept free and clean.
- Doors and drawers should only be open when items are being removed from or placed in the safety cabinet – keep the doors/drawers closed!
- Please observe the maximum amounts that can be stored, shelf loading, etc., as per Annex 1 and the markings on the door/drawer.
- The safety cabinets are not intended for the storage of aggressive substances (e.g. acids and alkalis, corrosive gases) as these could corrode the interior fittings and the exhaust air extraction system! Our special acid and alkali safety cabinets should be used for this. Attention! The storage of corrosive gases in gas cylinder cabinets can adversely affect the function of the supply and exhaust mechanism. This must be checked at regular intervals.
- Before storing substances (e.g. reagents) you must check that the cabinet surfaces are resistant to them.
- To comply with EN 14470-1 the collection trough must be able to contain a minimum volume of 10% of the contents of the containers stored above it, or at least 110% of the contents of the largest single container – do not store too large containers!
- See 7.4 perforated cover plates
- Pipe and cable penetrations shall be carried out professionally according to fire protection regulations and must not impair the function of the cabinet. After removing the cable and pipe penetrations, the holes of fire protection must be locked, so that the overall design receives the fire protection class. Failure to comply with the cabinet loses its approval. (see instructions inside of cabinet door)

5. Transport

storeLAB® safety cabinets must only be lifted from beneath by using a suitable means of transport; they must be transported upright. Tilting during lifting and transport must be avoided at all costs. Tilting could cause damage that could limit or destroy the functioning of the built-in safety devices. You must also avoid jolting the cabinet when unloading it..

The carrying surface of the hand truck, fork-lift truck or similar must have at least the same area as the base of the safety cabinet.

Type 90		
		
Type 30		
		

Please pick up the cabinet from the front with a hand pallet truck or something similar. This location surface has to be pushed through up to the rear panel (depth 540 mm)

6. Installation

Placing the safety cabinet

The safety cabinet must be installed, operated and maintained in such a way that employees are safeguarded from the risks of fires and explosions. See assembly instructions (inside cabinet door).

Relocating the safety cabinet

For preventing any transport damage on the doors during relocating of the safety cabinet security bolts have to be previously removed from the doors. After relocating and realignment of the safety cabinet bolts are going to be screwed in again (Please see assembly instructions point 3).

7. Interior fittings

Shelves

- Installed according to assembly instructions (inside cabinet door)

Pull-out shelves

- Already built into cabinets

Shelves

- Already built into cabinets

Please observe the maximum shelf load (see Annex 1)!

If the shelves are not installed in accordance with the assembly instructions then the maximum shelf load cannot be guaranteed.

Collection trough / Perforated cover plate

The collection troughs in the base of the storeLAB® safety cabinet or the pull-out shelves are intended to collect any liquid leakages inside the cabinet and are not additional support surfaces for stored containers (see EN 14470-1).

If a perforated cover plate is used for the collection trough / pull-out shelf then this can be used as additional support surfaces for stored containers.

8. Storage amounts / Storing substances together

Valid national regulations must be observed.

9. Ventilation

All safety cabinets have an air inlet and an air outlet connection. The outlet connection (suction side, connection diameters attached annex 1) is connected to the exhaust system. The exhaust air is led away through ducting to discharge into the atmosphere at a non-hazardous location. Inlet air is taken from inside the room. Please observe the valid national rules, regulations and laws concerning the stipulated rate of air exchange and exhaust air removal. The function of the ventilation should be checked at regular intervals. (see 13.2)

Safety cabinets as per EN 14470-1 can also be operated without connection to a ventilation system. (In Germany please observe the requirements of TRGS-510; elsewhere the relevant national rules, regulations and laws.)

Please remember the following points for operation without connection to a ventilation system:

- The cabinet must be grounded to avoid electrostatic charges.
- The area 2.5 m around the safety cabinet up to a height of 0.5 m above floor level is to be regarded as being a Zone 2 explosion- endangered area.
- In accordance with the relevant accident prevention regulations (BGV A8), this Zone 2 explosion- endangered area must be clearly and visibly marked with the following hazard symbol:

The locations of the air inlet and outlet connections are shown under Annex 1 (see drawings).

10. Grounding

In accordance with BGR 132 (D), all safety cabinets are equipped with a grounding connection on the cabinet roof as standard. If liquids are transferred within the cabinet then the interior of the cabinet must be additionally grounded.

11. Markings

a) Markings as per EN 14470-1 / 14470-2 / TRGS-510 / EN ISO 7010 / ISO 3864.

b) Itemization serial number: aaaabbbccddddd - eeee

a = item number b = year of construction c = calendar week d = current number e = key number

12. Decommissioning after a fire

Depending on the duration of the fire and the fire load, an explosive vapour-air mixture can be formed inside any safety cabinet. Before opening the storeLAB® safety cabinet all sources of ignition, naked flames, etc., must be removed from a distance of 10 m around the safety cabinet.

- Inform the fire department about the contents of the cabinet and the period mentioned below.
- Only open the cabinet approx. 24 h after the fire has been extinguished, or wait for a period that is at least six times as long as the duration of the fire.
- Check that the outer surface of the cabinet is cold. If it is still hand-warm then the waiting period must be extended.
- Avoid under all circumstances the use of cutting tools that produce sparks and burning tools.
- Ensure that suitable fire-fighting agents are at hand during the opening procedure.
- Wear suitable protective clothing!

13. Care and maintenance

The storeLAB® safety cabinet is a technical safety device according to the relevant German and European safety regulations (complies with EEC directive). To ensure its protective function in the event of a fire it must be checked at regular intervals. The German safety regulations stipulate that the safety devices must be checked by a specialist at regular intervals. The interval between safety checks is currently one year.

Daily checks

- Escaped liquids must be removed immediately and disposed of properly.
- Check the storeLAB® safety cabinet regularly for mechanical damage, dirt and corrosion.

Monthly checks

The monthly check covers the following points:

- | | |
|--------------------|---|
| 1. Door closing: | Open the door and check that it closes properly. |
| 2. Ventilation: | Test the efficiency of the ventilation system with a piece of wool, meter, etc. |
| 3. Sealing strips: | Check that the sealing strips in the door groove and the side(s) of the door(s) are correctly seated. If any damage is visible the sealing strips must be replaced immediately. |

14. Disposal

The storeLab® safety cabinet can be completely dismantled. Individual components, e.g. metal, insulation panels, door locks, etc. can be recycled separately. National and local waste disposal regulations must be observed. To preserve resources the safety cabinet components or the safety cabinet itself should not be disposed of as domestic waste.

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declare under our sole responsibility that the Product

storeLAB® *Sicherheitsschrank* Modell 1400 / Modell 1200 / Modell 900 /
Modell 600 / Modell SIS-FAS 900 / Modell SIS-FAS 1550 / Modell SIS-FAS 1550-1430
Modell SIS-Typ 30 600 / Modell SIS-Typ 30 600-GL / Modell SIS-Typ 30 1200 /
Modell Typ 30 1200-GL

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storeLAB® *Safety Storage-Cabinet* model 1400 / model 1200 / model 900 /

model 600 / model SIS-FAS 900 / model SIS-FAS 1550 / Model SIS-FAS 1550-1430

model SIS-Type 30 600 / model SIS-Type 30 600-GL / model Type 30 1200 /

model Type 30 1200-GL

Underbench cabinets model 1400 FT / model 1400 SF / model 1400 SFSF /
model 1400 FTSF /

model 1100 FT / model 1100 SF / model 1100 SFSF / model 1100 FTSF /

model 880 FT / model 880 SF / model 880 FT SF / model **880 SF** SF /

model 600 FT / model 600 SF

model 1400-500FT / model 1400-500 SF / model 1400-500 FT SF / model 1400-500SFSF

model 1100-500 FT / model 1100-500 SF / model 1100-500 FT SF /

model 1100-500 SF SF /

model 880-500 FT / model 880-500 SF / model 880-500 FT SF / model 880-500 SF SF /

model 600-500 FT / model 600-500 SF /

Safety cabinets for documents (DKS) model 1200 / model 1400

Chemical/ cabinets with safety box model **950-SiB** 30 / model 950 GL-SiB 30 /

model 950-SiB 60 / model 950 GL-SiB 60 / model 950-2FAS-SiB 30 / model 950-2FAS GL-
SiB 30 / model 950-2FAS-SiB 60 / model 950-2FAS GL-SiB 60

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to which this declaration relates is in conformity with the relevant provisions of the following standards
or other normative documents.

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Relevant EEC-
Directives

EN 292 Safety of machinery- Basic concepts, general principles
For design

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The operation company is responsible for installing the products in accordance to national regulations which may apply in respect to the installation, use and also for regular inspections.

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LA CONT Umwelttechnik GmbH